


DOUJET series

Pre-packed, closed system



DOUJET 60GL
DOUJET 80GL
DOUJET 40GF
DOUJET 60GF


Pouch system



DOUJET 40PL
DOUJET 40PF

DOUJET Accessories

DOUJET Injector



D1011

INJECTA Inc
7, Gongdan-ro 160beon-gil,
Gunpo-si, Gyeonggi-do 15845 Korea
Phone: +82-31-426-3636
Fax: +82-31-426-6336
www.injecta-hq.com

INJECTA



DOUJET
Powder-Liquid Bone Cement
containing Gentamicin sulfate

DOUJET consists of liquid and powder to be bone cement based on radiopaque polymethyl methacrylate and to contain aminoglycoside antibiotic gentamicin.


Why do we need a closed system ?

- **ACGIH** (American Conference of Governmental Industrial Hygienists)
Short Term Exposure Limit (STEL) : 100ppm TWA
15minute exposure to MMA not to be exceeded during a workday
- **OSHA** (Occupational Safety & Health Administration) *Hospital Investigations : Health Hazards*
"Mixing of MMA should be done in a closed system"




DOUJET L&F series

Setting time



DOUJET L series



DOUJET F series



DOUJET Pre-packed, closed system

Vacuum powder cartridge, Liquid ampoule, Vacuum tube, Liquid funnel, Injection nozzle

- Simplest vacuum mixing system
 - Clear mixing system
 - Perfect close mixing system
 - No cement lose
- Less infection - No air exposure during mixing
 - Economic all-in-one system